






# Bellman Neckloop

BE9159



# Technical Specifications

- Neckloop diameter:** 22 cm / 9 inches  
**Weight:** 60 g / 1.6 oz  
**Cable length:** 90 cm / 3 ft  
**Connectors:** 3.5 mm tele plug (break-away connector on loop cable)  
**Load impedance:**  $2 \times 5 \Omega$   
**Magnetic output:** 1500 mA/m at 15 cm / 6 inches distance  
(2 x 50 mW input signal)



**European Office**  
Bellman & Symfon Europe AB  
Södra Långebergsgatan 30  
421 32 Västra Frölunda  
Sweden  
Phone +46 31 68 28 20  
Fax +46 31 68 28 90  
E-mail [info@bellman.se](mailto:info@bellman.se)  
[www.bellman.com](http://www.bellman.com)

**Asia Pacific Office**  
Bellman & Symfon Asia Ltd.  
1115-16, Main Tower, GD Int'l Hotel  
339 Huanshi Dong Rd  
Guangzhou, China, 510095  
Phone +86 20 8318 2950  
Fax +86 20 8333 5528  
E-mail [info.cn@bellman.com](mailto:info.cn@bellman.com)  
[www.bellman.cn](http://www.bellman.cn)